

DrugCheck 3000

Saliva Drug Test Kit



Use with the DrugCheck app for convenient record keeping!



Obtain samples easily and safely

The DrugCheck 3000 consists of two components: an hygienic sample collector for obtaining a saliva sample, and a test cassette for the analysis. The test cassette contains the buffer liquid and a window with two test strips, which displays the control and test lines.

Three easy steps for testing

The drug test itself is performed in three easy steps: swab the saliva sample, shake the test kit, wait for a brief incubation period, and then start the test. A colour indicator on the swab disappears as soon as it has absorbed enough oral fluid for a test.

Unambiguous on-the-spot drug screening

Its compact pocket sized design makes the DrugCheck 3000 easy to transport. It can be made ready for use quickly and easily, and is suitable for flexible use applications such as roadside checks and testing in safety-relevant workplaces.

**Small Sample
Required**

**Verified to AS/NZS
4760:2019***

DrugCheck 3000 6 Panel with OXY

Tests for 6 Drugs

**Sample Sufficiency
Indicator**

Made in Germany

**Results In 3 - 5
Minutes**

**Record Test Results
with the Free
DrugCheck App**

**Bulk Buy Discounts
Available**



**Buy Online
Now!**



DrugCheck 3000

Saliva Drug Test Kit

Drugs Tested & Cutoff Levels:

Amphetamines (AMP)
- 50 ng/mL

Cocaine (COC)
- 50 ng/mL

Marijuana (THC)
- 15 ng/mL

Methamphetamine (MET)
- 50 ng/mL

Opiates (OPI)
- 50 ng/mL

+ Either

Oxycodone (OXY)
- 40 ng/mL

OR

Benzodiazepines (BZO)*
- 15 ng/mL

Please note: Only the OXY variant of the DrugCheck 3000 is verified & compliant to AS/NZS 4760:2019



Step 1: Collect saliva sample by placing collector in mouth and moving around cheeks and gums. Red indicator will disappear once sufficient sample has been collected.



Step 2: Insert the sample collector into the funnel shaped opening in the test cassette. Ensure there is no gap between the funnel opening and the sample collector




Step 3: Shake the test for at least 15 seconds until the second red indicator ring has disappeared. Then place the test kit on an even surface and wait 60 seconds



Step 4: Remove the security tab on the front of the cassette and forcefully press the test kit down onto the sample collector as far as it will go.



Step 5: Evaluate the results as they become available. Ensure the "C" line has formed on both test strips for a valid result. A not-negative (positive) result can be interpreted 5 minutes after starting step 4.



Step 6: Test results can be optionally saved and stored in the free DrugCheck app available for iPhone and Android devices.

(03) 9125 4520

drugtestwarehouse.com.au

sales@drugtestwarehouse.com.au



**Buy Online
Now!**